

SPLCoat™ Poly-D-Lysine – Cell Culture Slide

(#30304/30308)

1. Description

SPL Life Sciences provides Poly-D-Lysine coated cell culture slides with high quality. All Poly-D-Lysine coated slides are produced with optically clear and high quality polymers for microscopy.

2. Features

- ❖ Poly-D-Lysine (Source: Synthetic)
- ❖ No washing required before use
- ❖ Shelf life: 6 months at room temperature

3. Product Range

Cat.No.	Type	Growth Area (cm ²)	Working Volume (ml)	Packaging
30304	4well	2.13	0.50 – 1.30	6 / 12
30308	8well	0.98	0.20 – 0.60	6 / 12

4. Certification


SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Specifications & Technical Drawing

(1) 4well Cell Culture Slide (# 30304)
- Specifications

Dimensions			Material (Resin)	Packaging
Lid	Length:	24.59 mm	Polystyrene (PS)	▪ Pieces/Sleeve: 6 ▪ Pieces/Case: 12
	Width:	52.69 mm		
	Height:	7.10 mm		
Chamber	Length:	26.50 mm	Polystyrene (PS)	
	Width:	52.40 mm		
	Height:	11.30 mm		
Slide	Length:	26.00 mm	Glass	
	Width:	76.00 mm		
	Height:	1.00 mm		
Holder	Length:	27.58 mm	Polypropylene (PP)	
	Width:	56.00 mm		
	Height:	4.50 mm		
Well Dimension:		23.00 x 11.40 mm		
Volume per Well:		0.5 ~ 1.3 ml		


- Technical Drawing



(2) 8well Cell Culture Slide (# 30308)
 - Specifications

Dimensions			Material (Resin)	Packaging
Lid	Length:	24.59 mm	Polystyrene (PS)	▪ Pieces/Sleeve: 6 ▪ Pieces/Case: 12
	Width:	52.69 mm		
	Height:	7.10 mm		
Chamber	Length:	26.50 mm	Polystyrene (PS)	
	Width:	52.40 mm		
	Height:	11.30 mm		
Slide	Length:	26.00 mm	Glass	
	Width:	76.00 mm		
	Height:	1.00 mm		
Holder	Length:	27.58 mm	Polypropylene (PP)	
	Width:	56.00 mm		
	Height:	4.50 mm		
Well Dimension:		11.50 x 11.40 mm		
Volume per Well:		0.2 ~ 0.6 ml		

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
 To place an order, contact your local distributor or
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-CCS-PDL

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
 © SPL Life Sciences Co., Ltd., Version 7, November 20th, 2018

